

Product Name: SPT4H Rabbit Polyclonal Antibody
Catalog #: APRab18224



Summary

Production Name	SPT4H Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	SUPT4H1 SPT4H SUPT4H
Alternative Names	
Gene ID	6827.0
SwissProt ID	P63272.Synthesized peptide derived from human protein . at AA range: 40-120

Application

Dilution Ratio	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	12kD

Background

This gene encodes the small subunit of DRB (5,6-dichloro-1-beta-d-ribofuranosylbenzimidazole) sensitivity-inducing factor (DSIF) complex, which regulates mRNA processing and transcription elongation by RNA polymerase II. The encoded protein


Product Name: SPT4H Rabbit Polyclonal Antibody
Catalog #: APRab18224



is localized to the nucleus and interacts with the large subunit (SUPT5H) to form the DSIF complex. Related pseudogenes have been identified on chromosomes 2 and 12. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Nov 2012],function:Component of the DRB sensitivity-inducing factor complex (DSIF complex), which regulates mRNA processing and transcription elongation by RNA polymerase II. DSIF positively regulates mRNA capping by stimulating the mRNA guanylyltransferase activity of RNGTT/CAP1A. DSIF also acts cooperatively with the negative elongation factor complex (NELF complex) to enhance transcriptional pausing at sites proximal to the promoter. Transcriptional pausing may facilitate the assembly of an elongation competent RNA polymerase II complex. DSIF and NELF promote pausing by inhibition of the transcription elongation factor TFIIS/S-II. TFIIS/S-II binds to RNA polymerase II at transcription pause sites and stimulates the weak intrinsic nuclease activity of the enzyme. Cleavage of blocked transcripts by RNA polymerase II promotes the resumption of transcription from the new 3' terminus and may allow repeated attempts at transcription through natural pause sites. DSIF can also positively regulate transcriptional elongation and is required for the efficient activation of transcriptional elongation by the HIV-1 nuclear transcriptional activator, Tat. DSIF acts to suppress transcriptional pausing in transcripts derived from the HIV-1 LTR and blocks premature release of HIV-1 transcripts at terminator sequences.,similarity:Belongs to the SPT4 family.,subunit:Interacts with SUPT5H to form DSIF. DSIF interacts with the positive transcription elongation factor b complex (P-TEFb complex), which is composed of CDK9 and cyclin-T (CCNT1 or CCNT2). DSIF interacts with RNA polymerase II, and this interaction is reduced by phosphorylation of the C-terminal domain (CTD) of POLR2A by P-TEFb. DSIF also interacts with the NELF complex, which is composed of WHSC2/NELFA, COBRA1/NELFB, TH1L/NELFD and RDBP/NELFE, and this interaction occurs following prior binding of DSIF to RNA polymerase II. DSIF also interacts with HRMT1L2/PRMT1, HTATSF1/TATSF1, RNGTT/CAP1A, SKB1/PRMT5, SUPT6H, and can interact with PIN1.,tissue specificity:Widely expressed,.

Research Area

Image Data



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

Product Name: SPT4H Rabbit Polyclonal Antibody
Catalog #: APRab18224



Note

For research use only.