

Summary

Production Name	StARD10 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	STARD10 STARD10; SDCCAG28; CGI-52; PCTP-like protein; PCTP-L; Antigen NY-CO-28; START
Alternative Names	domain-containing protein 10; StARD10; Serologically defined colon cancer antigen 28; StAR-related lipid transfer protein 10
Gene ID	10809.0
SwissProt ID	Q9Y365.Synthesized peptide derived from the Internal region of human StARD10.

Application

Dilution Ratio	WB 1:500-2000 ELISA 2000-20000
Molecular Weight	33kD

Background


Product Name: StARD10 Rabbit Polyclonal Antibody
Catalog #: APRab18341



function:May play specific roles in sperm maturation or fertilization.,sequence caution:Various sequencing problems as well as a translation in a wrong frame.,similarity:Contains 1 START domain.,function:May play specific roles in sperm maturation or fertilization.,sequence caution:Various sequencing problems as well as a translation in a wrong frame.,similarity:Contains 1 START domain.,

Research Area

Image Data



Western Blot analysis of CoLo cells using StARD10 Polyclonal Antibody

Note

For research use only.