

Summary

Production Name	TCAL2 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	TCEAL2 My048 MY0876G0
Alternative Names	
Gene ID	140597.0
SwissProt ID	Q9H3H9.Synthesized peptide derived from part region of human protein

Application

Dilution Ratio	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	24kD

Background

This gene encodes a member of the transcription elongation factor A (SII)-like (TCEAL) gene family. Members of this family contain TFA domains and may function as nuclear phosphoproteins that modulate transcription in a promoter context-


Product Name: TCAL2 Rabbit Polyclonal Antibody
Catalog #: APRab18715



dependent manner. Multiple family members are located on the X chromosome. [provided by RefSeq, Jul 2008],function:May be involved in transcriptional regulation.,similarity:Belongs to the TFS-II family. TFA subfamily.,

Research Area

Image Data



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

Note

For research use only.