


Product Name: Recombinant Human FABP1 (N-6His)
Catalog #: PEH0614



Summary

Name	FABP1/L-FABP
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Fatty Acid-Binding Protein 1 is produced by our E.coli expression system and the target gene encoding Met1-Ile127 is expressed with a 6His tag at the N-terminus.
Accession #	P07148
Host	E.coli
Species	Human
Predicted Molecular Mass	16.37 KDa
Formulation	Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, 0.5M Argine, 50% Glycerol, 2mM EDTA, pH 7.4.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image



Background

Alternative Names	Fatty Acid-Binding Protein Liver; Fatty Acid-Binding Protein 1; Liver-Type Fatty Acid-Binding Protein; L-FABP; FABP1; FABPL
Background	Human Fatty Acid-Binding Protein 1 (FABP1) is a cytoplasm protein, which belongs to the calycin superfamily and Fatty-acid binding protein (FABP) family. Fatty acid

Product Name: Recombinant Human FABP1 (N-6His)
Catalog #: PEH0614



binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids. FABP1 forms a beta-barrel structure that accommodates hydrophobic ligands in its interior. FABP1 can bind free fatty acids and their coenzyme A derivatives, bilirubin, and some other small molecules in the cytoplasm, so it can be involved in intracellular lipid transport.

Note

For Research Use Only , Not for Diagnostic Use.