


Product Name: Recombinant Human MIF(N-6His)
Catalog #: PEH1164



Summary

Name	MIF/Macrophage migration inhibitory factor
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Macrophage Migration Inhibitory Factor is produced by our E.coli expression system and the target gene encoding Met1-Ala115 is expressed with a 6His tag at the N-terminus.
Accession #	P14174
Host	E.coli
Species	Human
Predicted Molecular Mass	14.6 KDa
Formulation	Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, 50% Glycerol, pH7.4.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	0.00.0

SDS-PAGE image



Background

Alternative Names	Macrophage migration inhibitory factor; MIF; MMIF; Glycosylation-inhibiting factor; GLIF; L-dopachrome tautomerase; Phenylpyruvate tautomerase
Background	Human MIF is a 12.5 kDa, 115 amino acid (aa) nonglycosylated polypeptide that is synthesized without a signal sequence. Secretion occurs nonclassically via an

Product Name: Recombinant Human MIF(N-6His)
Catalog #: PEH1164



ABCA1 transporter. Pro-inflammatory cytokine. Involved in the innate immune response to bacterial pathogens. The expression of MIF at sites of inflammation suggests a role as mediator in regulating the function of macrophages in host defense. Counteracts the anti-inflammatory activity of glucocorticoids. Has phenylpyruvate tautomerase and dopachrome tautomerase activity (in vitro), but the physiological substrate is not known. It is not clear whether the tautomerase activity has any physiological relevance, and whether it is important for cytokine activity.

Note

For Research Use Only , Not for Diagnostic Use.