


Summary

Name	peroxiredoxin 6/PRDX6/AOP2
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Peroxiredoxin-6 is produced by our E.coli expression system and the target gene encoding Met1-Pro224 is expressed with a 6His tag at the N-terminus.
Accession #	P30041
Host	E.coli
Species	Human
Predicted Molecular Mass	27.2 KDa
Formulation	Supplied as a 0.2 μm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image



Background

Alternative Names	Peroxiredoxin-6; 1-Cys Peroxiredoxin; 1-Cys PRX; 24 kDa Protein; Acidic Calcium-Independent Phospholipase A2; aiPLA2; Antioxidant Protein 2; Liver 2D PAGE Spot 40; Non-Selenium Glutathione Peroxidase; NSGPx; Red Blood Cells Page Spot 12; PRDX6; AOP2; KIAA0106
--------------------------	---

Product Name: Recombinant Human PRDX6 (N-6His)
Catalog #: PEH1308



Background

Peroxiredoxin 6 belongs to the thiol-specific antioxidant protein family, which is a bifunctional enzyme with two distinct active sites. Peroxiredoxin 6 is involved in redox regulation of the cell and can reduce Hydrogen peroxide and short chain organic, fatty acid, and phospholipid hydroperoxides. Peroxiredoxin 6 may regulate phospholipid turnover and protect against oxidative injury. Peroxiredoxin 6 eases the oxidative stress and TGF-beta-induced abnormalities of human trabecular meshwork cells, it is necessary for blood vessel integrity in injured skin.

Note

For Research Use Only , Not for Diagnostic Use.