


Product Name: Recombinant Human TXN (N-6His)
Catalog #: PEH1632



Summary

Name	Thioredoxin/TXN/Trx
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Thioredoxin is produced by our E.coli expression system and the target gene encoding Met1-Val105 is expressed with a 6His tag at the N-terminus.
Accession #	P10599
Host	E.coli
Species	Human
Predicted Molecular Mass	13.9 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 1mM EDTA, 2mM DTT, pH 7.2.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Background

Product Name: Recombinant Human TXN (N-6His)
Catalog #: PEH1632



Alternative Names

Thioredoxin; Trx; ATL-Derived Factor; ADF; Surface-Associated Sulphydryl Protein; SASP; TXN; TRDX; TRX; TRX1

Background

Thioredoxin (TXN) is a member of the Thioredoxin family. Thioredoxin exists as a disulfide-linked homodimer and contains one Thioredoxin domain. Thioredoxin is up-regulated by ionizing radiation. Thioredoxin participates in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyzes dithiol-disulfide exchange reactions. Thioredoxin also plays a role in the reversible S-nitrosylation of cysteine residues in target proteins, and thereby contributes to the response to intracellular nitric oxide.

Note

For Research Use Only , Not for Diagnostic Use.