


Summary

Name	CD158d/KIR2DL4
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Killer Cell immunoglobulin-like Receptor 2DL4 is produced by our Mammalian expression system and the target gene encoding Trp22-His242 is expressed with a 6His tag at the C-terminus.
Accession #	ADY38409.1
Host	Human Cells
Species	Human
Predicted Molecular Mass	25.34 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Product Name: Recombinant Human KIR2DL4 (C-6His)
Catalog #: PHH0296



Background

Alternative Names

Killer Cell Immunoglobulin-Like Receptor 2DL4; CD158 Antigen-Like Family Member D; G9P; Killer Cell Inhibitory Receptor 103AS; KIR-103AS; MHC Class I NK Cell Receptor KIR103AS; CD158d; KIR2DL4; CD158D; KIR103AS

Background

Killer cell immunoglobulin-like receptor 2DL4(KIR2DL4) is a Single-pass type I membrane protein and contains 2 Ig-like C2-type (immunoglobulin-like) domains. It belongs to the immunoglobulin superfamily. KIR2DL4 is expressed in all NK cells and some T cells. KIR2DL4 activates the cytotoxicity of NK cells, despite the presence of an immunoreceptor tyrosine-based inhibition motif (ITIM) in its cytoplasmic tail. The ITIM was not necessary for activation of lysis by KIR2DL4. The activation signal of KIR2DL4 was sensitive to inhibition by another ITIM-containing receptor. The activation-deficient mutant of KIR2DL4 inhibited the signal delivered by the activating receptor CD16.

Note

For Research Use Only , Not for Diagnostic Use.