


Summary

Name	Secretogranin-3/SCG3
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Secretogranin-3 is produced by our Mammalian expression system and the target gene encoding Phe20-Leu468 is expressed with a 6His tag at the C-terminus.
Accession #	AAH14539.1
Host	Human Cells
Species	Human
Predicted Molecular Mass	52 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Background

Product Name: Recombinant Human SCG3 (C-6His)
Catalog #: PHH1476



Alternative Names Secretogranin-3;Secretogranin III;SgIII;

Background Secretogranin-3(SCG3) is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. SCG3 is expressed in brain, heart, kidney, liver and skeletal muscle. It may serve as precursors for biologically active peptides. Some granins have been shown to function as helper proteins in sorting and proteolytic processing of prohormones; however, the function of this protein is unknown.

Note

For Research Use Only , Not for Diagnostic Use.