


Summary

Name	Thiamin pyrophosphokinase 1/TPK1
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Thiamin Pyrophosphokinase 1 is produced by our Mammalian expression system and the target gene encoding Met1-Ser243 is expressed with a 6His tag at the C-terminus.
Accession #	AAH68460.1
Host	Human Cells
Species	Human
Predicted Molecular Mass	28.3 KDa
Formulation	Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image



Background

Alternative Names	Thiamin pyrophosphokinase 1; hTPK1; Placental protein 20; PP20; Thiamine pyrophosphokinase 1; TPK1
Background	Thiamin pyrophosphokinase 1 (TPK1), belongs to the thiamine pyrophosphokinase family. TPK1 exists as a homodimer, and is detected in heart, kidney, testis, small

Product Name: Recombinant Human TPK1 (C-6His)
Catalog #: PHH1628



intestine and peripheral blood leukocytes, and at very low levels in a variety of tissues. TPK1 can catalyze the phosphorylation of thiamine to thiamine pyrophosphate. It can also catalyze the phosphorylation of pyriothiamine to pyriothiamine pyrophosphate.

Note

For Research Use Only , Not for Diagnostic Use.