


Summary

Name	TrkB/NTRK2/Neurotrophic Tyrosine Kinase Receptor Type 2
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/ μ g as determined by LAL test.
Construction	Recombinant Human Neurotrophic Tyrosine Kinase Receptor Type 2 is produced by our Mammalian expression system and the target gene encoding Cys32-His430 is expressed with a 6His tag at the C-terminus.
Accession #	Q16620
Host	Human Cells
Species	Human
Predicted Molecular Mass	45.25 KDa
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at \leq -70°C, stable for 6 months after receipt. Store at \leq -70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μ g/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Background

Product Name: Recombinant Human TrkB (C-6His)
Catalog #: PHH1728



Alternative Names

BDNF/NT-3 Growth Factors Receptor; GP145-TrkB; Trk-B; Neurotrophic Tyrosine Kinase Receptor Type 2; TrkB Tyrosine Kinase; Tropomyosin-Related Kinase B; NTRK2; TRKB

Background

The TRK Family of Tyrosine Kinase Receptor consists of 3 members: TrkA, TrkB and TrkC. The three TRK family proteins have different ligand specificities. They connect to different neurotrophins, including NGF, BDNF, NT-3/NT-4/5. TRKA binds NGF, TRKB binds BDNF and NT-3, TRKC binds NT-4/5. At the protein sequence level, human and rat TRKB have greater than 90% sequence identity and the proteins exhibit cross-species activity. TRKB is primarily expressed in the nervous system and it also expression in a wide variety of tissues with low levels.

Note

For Research Use Only , Not for Diagnostic Use.