


Summary

Name	FGL-1/Fibrinogen-like protein 1
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Fibrinogen-like Protein 1 is produced by our Mammalian expression system and the target gene encoding Leu23-Ile312 is expressed with a 6His tag at the C-terminus.
Accession #	Q08830
Host	Human Cells
Species	Human
Predicted Molecular Mass	35 KDa
Formulation	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 300mM NaCl, 10% Glycerol, 2mM EDTA, pH 8.0.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image



Background

Alternative Names Fibrinogen-like protein 1; FGL1; HP-041; Hepassocin; HFREP-1; LFIRE-1

Background Fibrinogen-like protein 1(FGL1) is also known as HP-041, Hepassocin, HFREP-1, LFIRE-1, is a liver-specific secreted protein belonging to the fibronogen

Product Name: Recombinant Human FGL1 (C-6His)
Catalog #: PHH2064



superfamily, whose members share a fibrinogen domain at their C-termini. It is secreted by the liver and functions as a mitogen for hepatocytes. Hepassocin may play a role in the development of hepatocellular carcinomas. Hepassocin is a disulfide-linked homodimeric protein with a C-terminal fibrinogen domain. It is reported that it is a major immune inhibitory ligand of LAG-3 .

Note

For Research Use Only , Not for Diagnostic Use.