


Summary

Name	GAS6
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Growth Arrest-specific Protein 6 is produced by our Mammalian expression system and the target gene encoding Ala31-Ala678 is expressed with a 6His tag at the C-terminus.
Accession #	Q14393-2
Host	Human Cells
Species	Human
Predicted Molecular Mass	72.7 KDa
Formulation	Supplied as a 0.2 μm filtered solution of PBS, 10% Glycerol, pH 7.4.
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	

SDS-PAGE image



Background

Alternative Names	AXLLG; AXLLGAXL stimulatory factor; AXSFAXL receptor tyrosine kinase ligand; Gas6; GAS-6; growth arrest-specific 6; growth arrest-specific protein 6
Background	GAS6 (Growth arrest-specific protein 6) is also known as AXL receptor tyrosine

Product Name: Recombinant Human GAS6 (C-6His)
Catalog #: PHH2163



kinase ligand, AXLLG, is a multimodular protein that is up-regulated by a wide variety of cell types in response to growth arrest. Gas6 binds and induces signaling through the receptor tyrosine kinases Axl, Dtk, and Mer whose signaling is implicated in cell growth and survival, cell adhesion and cell migration. GAS6/AXL signaling plays a role in various processes such as endothelial cell survival during acidification by preventing apoptosis, optimal cytokine signaling during human natural killer cell development, hepatic regeneration, gonadotropin-releasing hormone neuron survival and migration, platelet activation, or regulation of thrombotic responses.

Note

For Research Use Only , Not for Diagnostic Use.