


Summary

Name	TACI/TNFRSF13B/CD267
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Transmembrane Activator and CAML Interactor is produced by our Mammalian expression system and the target gene encoding Ser2/xadThr166 is expressed with a human IgG1 Fc tag at the C-terminus.
Accession #	O14836
Host	Human Cells
Species	Human
Predicted Molecular Mass	45.6 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Background

Product Name: Recombinant Human TACI (C-Fc)
Catalog #: PHH2270



Alternative Names

Tumor necrosis factor receptor superfamily member 13B; Transmembrane activator and CAML interactor; CD267; Tnfrsf13b

Background

Tumor necrosis factor receptor superfamily member 13B, also known as TNFRSF13B or more commonly as TACI (transmembrane activator and CAML interactor), is a transmembrane receptor protein found predominantly on the surface of B cells, which are an important part of the immune system.

Note

For Research Use Only , Not for Diagnostic Use.