


Product Name: Recombinant Human Laminin521 E8 Protein
Catalog #: PHH2461



Summary

Name	Laminin521 E8
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Construction	Recombinant Human Laminin521 is produced by our Mammalian expression system and the target gene encoding Ala2534-Ala3322&Leu1573-Gln1798&Asn1364-Pro1609 is expressed with a His, Flag and HA tag at N terminus respectively.
Accession #	O15230 & P55268 & P11047
Host	Human cells
Species	Human
Predicted Molecular Mass	86.7&25.8&28.5 Kda
Formulation	0.0
Shipping	0.0
Stability&Storage	0.0
Reconstitution	0.00.0

SDS-PAGE image



Background

Product Name: Recombinant Human Laminin521 E8 Protein
Catalog #: PHH2461



Alternative Names 0.0

Background 0.0

Note

For Research Use Only , Not for Diagnostic Use.