

**Product Name:** Monoclonal Mouse Anti-Human IgG(Fab specific)

**Catalog #:** AMS0002

---

## Summary

<b>Production Name</b>	Monoclonal Mouse Anti-Human IgG(Fab specific)
<b>Clonality</b>	Monoclonal
<b>Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Application</b>	ELISA,WB
<b>Form</b>	Liquid
<b>Concentration</b>	5mg/mL
<b>Immunogen</b>	IgG
<b>Specificity</b>	IgG(Fab specific)
<b>Buffer</b>	0.01M PBS、 PH:7.2、 N preservative
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodification
<b>Purification</b>	Affinity purification
<b>Purity</b>	Calculated as antibodies, no less than 90%.
<b>Storage</b>	Please store at -20°C away from light.

## Application

ELISA: 1:20000-1:100000, WB: 1:20000-50000. The clone number of this monoclonal antibody is HG-12, the monoclonal antibody is mouse IgG1 subclass, and the light chain is Kappa. It is suitable for labeling and ELISA experiments.

## Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.