

Product Name: Polyclonal Rabbit Anti-Mouse IgG(H+L)
Catalog #: APS0006



Summary

| | |
|------------------------|---|
| Production Name | Polyclonal Rabbit Anti-Mouse IgG(H+L) |
| Clonality | Polyclonal |
| Host | Rabbit |
| Reactivity | Mouse |
| Application | ELISA,ICC/IF,Flow Cyt,WB |
| Form | Liquid |
| Concentration | 5mg/mL |
| Immunogen | IgG |
| Specificity | IgG(H+L) |
| Buffer | 0.01M PBS、 PH:7.2、 N preservative |
| Conjugation | Unconjugated |
| Modification | Unmodification |
| Purification | Affinity purification |
| Purity | Calculated as antibodies, no less than 90%. |
| Storage | Please store at -20°C away from light. |

Application

ELISA: 1:20000-1:100000; ICC/IF: 1:200-1:2000; Flow Cyt: 1:1000-4000; WB: 1:20000-50000; This data is a recommended range, please refer to the actual value according to the specific situation. This product can also be well applied to conventional in vitro immunological experiments such as colloidal gold.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.