

Product Name: Polyclonal Goat Anti-Rabbit IgG (H+L)
(Non cross with Mouse IgG)
Catalog #: APS0015



Summary

Production Name	Polyclonal Goat Anti-Rabbit IgG (H+L)(Non cross with Mouse IgG)
Clonality	Polyclonal
Host	Goat
Reactivity	Rabbit
Application	ELISA,ICC/IF,Flow Cyt,WB
Form	Liquid
Concentration	5mg/mL
Immunogen	IgG
Specificity	IgG (H+L)(Non cross with Human IgG)
Buffer	0.01M PBS、 PH:7.2、 N preservative
Conjugation	Unconjugated
Modification	Unmodification
Purification	Affinity purification
Purity	Calculated as antibodies, no less than 90%.
Storage	Please store at -20°C away from light.

Application

ELISA: 1:20000-1:100000; ICC/IF: 1:200-1:2000; Flow Cyt: 1:1000-4000; WB: 1:20000-50000; This data is a recommended range, please refer to the actual value according to the specific situation. This product can also be well applied to conventional in vitro immunological experiments such as colloidal gold.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.