

Product Name: HRP-Polyclonal Goat Anti-Human IgM(μ -chain specific)

Catalog #: APS0027

Summary

Production Name	HRP-Polyclonal Goat Anti-Human IgM(μ -chain specific)
Clonality	Polyclonal
Host	Goat
Reactivity	Human
Application	ELISA,WB,IHC
Form	Liquid
Concentration	1mg/mL
Immunogen	IgM
Specificity	IgM(μ -chain specific)
Buffer	50% Glycerol, 0.01M PBS, Proclin-300
Conjugation	HRP
Modification	Unmodification
Purification	Affinity purification
Purity	Calculated as antibodies, no less than 90%.
Storage	Please store at -20°C away from light.

Application

ELISA: 1:4000-20000, WB: 1:4000-12000, IHC: 1:500-2500. This data is a recommended range, please refer to the actual value according to the specific situation. The recommended diluent is 0.01M PBS (PH: 7.0) containing 0.2% BSA and 0.05% Tween-20. Other suitable diluents can also be used. In principle, the prepared working solution should be used up once. If it needs to be stored for a long time and used multiple times after dilution, please contact us to purchase stabilizers or ready-to-use diluents.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.