

**Product Name:** HRP-Polyclonal Goat Anti-Mouse IgG (H+L)(Non cross with Rabbit IgG)

**Catalog #:** APS0029



## Summary

<b>Production Name</b>	HRP-Polyclonal Goat Anti-Mouse IgG (H+L)(Non cross with Rabbit IgG)
<b>Clonality</b>	Polyclonal
<b>Host</b>	Goat
<b>Reactivity</b>	Mouse
<b>Application</b>	ELISA,WB,IHC
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/mL
<b>Immunogen</b>	IgG
<b>Specificity</b>	IgG (H+L)(Non cross with Rabbit IgG)
<b>Buffer</b>	50% Glycerol, 0.01M PBS, Proclin-300
<b>Conjugation</b>	HRP
<b>Modification</b>	Unmodification
<b>Purification</b>	Affinity purification
<b>Purity</b>	Calculated as antibodies, no less than 90%.
<b>Storage</b>	Please store at -20°C away from light.

## Application

ELISA: 1:4000-20000, WB: 1:4000-12000, IHC: 1:500-2500. This data is a recommended range, please refer to the actual value according to the specific situation. The recommended diluent is 0.01M PBS (PH: 7.0) containing 0.2% BSA and 0.05% Tween-20. Other suitable diluents can also be used. In principle, the prepared working solution should be used up once. If it needs to be stored for a long time and used multiple times after dilution, please contact us to purchase stabilizers or ready-to-use diluents.

## Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.