

Product Name: FITC-Polyclonal Goat Anti-Mouse IgG (H+L)(Non cross with Rabbit IgG)

Catalog #: APS0044



Summary

Production Name	FITC-Polyclonal Goat Anti-Mouse IgG (H+L)(Non cross with Rabbit IgG)
Clonality	Polyclonal
Host	Goat
Reactivity	Mouse
Application	Flow-Cyt/IF,IHC
Form	Liquid
Concentration	1mg/mL
Immunogen	IgG
Specificity	IgG (H+L)(Non cross with Rabbit IgG)
Buffer	50% Glycerol, 0.01M PBS, 0.5%GSA, N preservative
Conjugation	FITC
Modification	Unmodification
Purification	Affinity purification
Purity	Calculated as antibodies, no less than 90%.
Storage	Please store at -20°C away from light.

Application

Flow-Cyt/IF:1:150-1500, IHC: 1:500-2500, This data is a recommended range, please refer to the actual value according to the specific situation. This product also can be used in Elisa experiment.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.