

Product Name: BIOTIN-Polyclonal Rabbit Anti-Human IgM(μ -chain specific)
Catalog #: APS0052



Summary

Production Name	BIOTIN-Polyclonal Rabbit Anti-Human IgM(μ -chain specific)
Clonality	Polyclonal
Host	Rabbit
Reactivity	Human
Application	Flow-Cyt/IF/IP,IHC,ELISA,WB
Form	Liquid
Concentration	1mg/mL
Immunogen	IgM
Specificity	IgM(μ -chain specific)
Buffer	50% Glycerol, 0.01M PBS, 0.5%RSA, N preservative
Conjugation	BIOTIN
Modification	Unmodification
Purification	Affinity purification
Purity	Calculated as antibodies, no less than 90%.
Storage	Please store at -20°C away from light.

Application

Flow-Cyt/IF/IP:1:200-2000, IHC: 1:500-2500, ELISA/WB: 1: 2000-20000. This data is a recommended range, please refer to the actual value according to the specific situation.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.