

**Product Name:** BIOTIN-Polyclonal Goat Anti-Mouse IgG (H+L)(Non cross with Rabbit IgG)  
**Catalog #:** APS0058

---

## Summary

<b>Production Name</b>	BIOTIN-Polyclonal Goat Anti-Mouse IgG (H+L)(Non cross with Rabbit IgG)
<b>Clonality</b>	Polyclonal
<b>Host</b>	Goat
<b>Reactivity</b>	Mouse
<b>Application</b>	Flow-Cyt/IF/IP,IHC,ELISA,WB
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/mL
<b>Immunogen</b>	IgG
<b>Specificity</b>	IgG (H+L)(Non cross with Rabbit IgG)
<b>Buffer</b>	50% Glycerol, 0.01M PBS, 0.5%GSA, N preservative
<b>Conjugation</b>	BIOTIN
<b>Modification</b>	Unmodification
<b>Purification</b>	Affinity purification
<b>Purity</b>	Calculated as antibodies, no less than 90%.
<b>Storage</b>	Please store at -20°C away from light.

## Application

Flow-Cyt/IF/IP:1:200-2000, IHC: 1:500-2500, ELISA/WB: 1: 2000-20000. This data is a recommended range, please refer to the actual value according to the specific situation.

## Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.