

Product Name: BIOTIN-Polyclonal Goat Anti-Rabbit IgG (H+L)(Non cross with Mouse IgG)
Catalog #: APS0059

Summary

| | |
|------------------------|--|
| Production Name | BIOTIN-Polyclonal Goat Anti-Rabbit IgG (H+L)(Non cross with Mouse IgG) |
| Clonality | Polyclonal |
| Host | Goat |
| Reactivity | Rabbit |
| Application | Flow-Cyt/IF/IP,IHC,ELISA,WB |
| Form | Liquid |
| Concentration | 1mg/mL |
| Immunogen | IgG |
| Specificity | IgG (H+L)(Non cross with Mouse IgG) |
| Buffer | 50% Glycerol, 0.01M PBS, 0.5%GSA, N preservative |
| Conjugation | BIOTIN |
| Modification | Unmodification |
| Purification | Affinity purification |
| Purity | Calculated as antibodies, no less than 90%. |
| Storage | Please store at -20°C away from light. |

Application

Flow-Cyt/IF/IP:1:200-2000, IHC: 1:500-2500, ELISA/WB: 1: 2000-20000. This data is a recommended range, please refer to the actual value according to the specific situation.

Note

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.