

Product Name: Mouse anti-Crucian carp IgM Heavy chain Monoclonal Antibody
Catalog #: AMM40008

Summary

Production Name	Mouse anti-Crucian carp IgM Heavy chain Monoclonal Antibody
Host	Mouse
Application	ELISA, WB, FC
Reactivity	Crucian carp

Performance

Conjugation	Unconjugation
Modification	Unmodification
Isotype	IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at -20°C. Stable for one year after shipment.
Buffer	1×PBS, 50% glycerin, 3 mg/mL BSA, pH 7.4.


Application

Recommended concentration	ELISA/WB:1:1000 ;FC:1:500
----------------------------------	---------------------------


Product Name: Mouse anti-Crucian carp IgM Heavy chain Monoclonal Antibody
Catalog #: AMM40008

Selected Validation Data


WB



ELISA




FC



IgM

Elution volume 6.5 7.5 8.5 9.5 10.5 11.5 12.5 13.5 14.5 (mL)



Note

For research use only.