

Product Name: DiMethyl-Histone H3 (Lys9) (10B2)
Mouse Monoclonal Antibody
Catalog #: AMM03582

Summary

Production Name	DiMethyl-Histone H3 (Lys9) (10B2) Mouse Monoclonal Antibody
Description	Primary antibody
Host	Mouse
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Methylated
Isotype	IgG1
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	H3C1
Alternative Names	H3K9me2; H3 histone; HIST1H3A; Histone cluster 1; H3a
Gene ID	8350
SwissProt ID	P68431

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 15 kDa; Observed MW: 15 kDa

Background


Product Name: DiMethyl-Histone H3 (Lys9) (10B2)
Mouse Monoclonal Antibody
Catalog #: AMM03582

H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

Research Area


Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of DiMethyl-Histone H3 in HeLa, rat Heart, Raw264.7 lysates using DiMethyl-Histone H3 (Lys9) (10B2)

antibody



Western blot analysis of DiMethyl-Histone H3 (Lys9) (10B2) in MCF-7, C2C12, C6, 3T3, HeLa lysates using DiMethyl-Histone H3 (Lys9) (10B2) antibody.

Note

For research use only.