

Summary

| | |
|------------------------|--|
| Production Name | ARFRP1 Rabbit Monoclonal Antibody |
| Description | Recombinant Rabbit Monoclonal antibody |
| Host | Rabbit |
| Application | WB,IP |
| Reactivity | Rat |

Performance

| | |
|---------------------|--|
| Conjugation | Unconjugated |
| Modification | Unmodified |
| Isotype | IgG |
| Clonality | Monoclonal Antibody |
| Form | Liquid |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Buffer | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA |
| Purification | Affinity Purified |

Immunogen

| | |
|--------------------------|------------------|
| Gene Name | ARFRP1 |
| Alternative Names | ARP; Arp1; ARL18 |
| Gene ID | 10139 |
| SwissProt ID | Q13795 |

Application

| | |
|-------------------------|--|
| Dilution Ratio | WB: 1/500-1/1000 IP: 1/20 |
| Molecular Weight | Calculated MW: 23 kDa; Observed MW: 23 kDa |

Background

Product Name: ARFRP1 Rabbit Monoclonal Antibody
Catalog #: AMRe01438




Trans-Golgi-associated GTPase that regulates protein sorting.

Research Area

Signal Transduction

Image Data



Western blot analysis of ARFRP1/ARP in C6 lysates using ARFRP1 antibody.

Note

For research use only.