

Product Name: Cytokeratin 4 Rabbit Monoclonal Antibody
Catalog #: AMRe01486

Summary

Production Name	Cytokeratin 4 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	KRT4
Alternative Names	Keratin4; Cytokeratin 4; Keratin-4; CK-4; CK4; KRT14; CYK4
Gene ID	3851
SwissProt ID	P19013

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 56 kDa

Product Name: Cytokeratin 4 Rabbit Monoclonal Antibody
Catalog #: AMRe01486


Background

Cytokeratins are a subfamily of intermediate filament keratins that are characterized by a remarkable biochemical diversity, which is represented in human epithelial tissues by at least 20 different polypeptides.

Research Area

Signal Transduction

Image Data



Western blot analysis of Cytokeratin 4 in C6, HeLa lysates using Cytokeratin 4 antibody.

Note

For research use only.