

Summary

Production Name	DKK3 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	DKK3
Alternative Names	DKK3; REIC; Dickkopf-related protein 3; Dickkopf-3; Dkk-3; hDkk-3
Gene ID	27122
SwissProt ID	Q9UBP4

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 70 kDa

Background

Product Name: DKK3 Rabbit Monoclonal Antibody
Catalog #: AMRe01914




Antagonizes canonical Wnt signaling by inhibiting LRP5/6 interaction with Wnt and by forming a ternary complex with the transmembrane protein KREMEN that promotes internalization of LRP5/6.

Research Area

Signal Transduction

Image Data



Western blot analysis of Dkk3 in rat Brain lysates using DKK3 antibody.

Note

For research use only.