
Summary

Production Name	DPS1 Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	PDSS1
Alternative Names	Trans-prenyltransferase 1; TPT 1; Solanesyl-diphosphate synthase subunit 1; Decaprenyl-diphosphate synthase subunit 1
Gene ID	23590
SwissProt ID	Q5T2R2

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW: 46 kDa; Observed MW: 46 kDa

Product Name: DPS1 Rabbit Monoclonal Antibody
Catalog #: AMRe01925




Background

Decaprenyl-diphosphate synthase subunit 1

Research Area

Signal Transduction

Image Data



Western blot analysis of DPS1 in K562 lysates using DPS1 antibody.

Note

For research use only.