

Product Name: HNF 4 alpha Rabbit Monoclonal Antibody
Catalog #: AMRe02102

Summary

Production Name	HNF 4 alpha Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP,ChIP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	Hnf4A
Alternative Names	Hnf4; HNF-4; MODY1; Nr2a1; Tcf14; TCF-14; Hnf4alpha
Gene ID	15378.0
SwissProt ID	P49698

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20 ChIP: 1/20
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 53 kDa

Product Name: HNF 4 alpha Rabbit Monoclonal Antibody
Catalog #: AMRe02102


Background

Transcriptionally controlled transcription factor. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha. May be essential for development of the liver, kidney and inTine.MiscellaneousBinds fatty acids.


Research Area

Cardiovascular

Image Data



Western blot analysis of HNF4alpha in mouse liver lysates using HNF 4 alpha antibody.



Western blot analysis of HNF4alpha in K562, C6, HeLa lysates using HNF4alpha antibody.

Note

For research use only.