

---

## Summary

<b>Production Name</b>	RanGAP1 Rabbit Monoclonal Antibody
<b>Description</b>	Recombinant Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC-F,IHC-P,ICC/IF,IP
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	RANGAP1
<b>Alternative Names</b>	SD; Fug1; RANGAP
<b>Gene ID</b>	5905
<b>SwissProt ID</b>	P46060

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20
<b>Molecular Weight</b>	Calculated MW: 64 kDa; Observed MW: 68,82 kDa

## Background

---

**Product Name: RanGAP1 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe03179**




GTPase activator for the nuclear Ras-related regulatory protein Ran, converting it to the putatively inactive GDP-bound state.


## Research Area

Cell Biology


## Image Data



Western blot analysis of RanGAP1 in 3T3 lysates using RanGAP1 antibody.




Western blot analysis of RanGAP1 in HeLa, A549, HL-60, U2OS, C6 lysates using RanGAP1 antibody.



Immunocytochemistry analysis of RanGAP1 (green) in HeLa using RanGAP1 antibody, and DAPI (blue)

**Product Name: RanGAP1 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe03179**



Immunohistochemistry analysis of paraffin-embedded Human tonsil using RanGAP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

**Note**

For research use only.