

Product Name: AGTR2 Receptor Rabbit Monoclonal Antibody
Catalog #: AMRe04001



Summary

Production Name	AGTR2 Receptor Rabbit Monoclonal Antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purified

Immunogen

Gene Name	AGTR2
Alternative Names	AGTR2; Angiotensin ii receptor type 2; Angiotensin II type-2 receptor; AT2 receptor; AT2; At2r; MRX88; Angiotensin ii at2; Angiotensin receptor 2; Type-2 angiotensin II receptor; ATGR2
Gene ID	186
SwissProt ID	P50052

Application

Dilution Ratio	WB: 1/500-1/1000 IP: 1/20
Molecular Weight	Calculated MW:41 kDa;Observed MW: 41 kDa

Product Name: AGTR2 Receptor Rabbit Monoclonal Antibody
Catalog #: AMRe04001


Background

Receptor for angiotensin II. Cooperates with MTUS1 to inhibit ERK2 activation and cell proliferation.

Research Area

Signal Transduction

Image Data



Western blot analysis of AGTR2 Receptor in HepG2, 293T, HeLa lysates using AGTR2 Receptor antibody.

Note

For research use only.