

Product Name: Cystatin B (7U18) Rabbit Monoclonal Antibody
Catalog #: AMRe09692



Summary

Production Name	Cystatin B (7U18) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	CSTB
Alternative Names	CHROW21; CST6; CSTB; Cystatin B (stefin B); CYTB EPM1; EPM1A; PME; STFB; ULD;
Gene ID	1476.0
SwissProt ID	P04080.A synthetic peptide of human Cystatin-B

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	11kDa

Background


Product Name: Cystatin B (7U18) Rabbit Monoclonal Antibody
Catalog #: AMRe09692



This is an intracellular thiol proteinase inhibitor. Tightly binding reversible inhibitor of cathepsins L, H and B. This is an intracellular thiol proteinase inhibitor. Tightly binding reversible inhibitor of cathepsins L, H and B.

Research Area

Image Data



Western blot detection of Cystatin-B in HeLa,A549,HL-60 using Cystatin-B antibody(1:1000 diluted)

Note

For research use only.