

Product Name: HSPA12A (5U3) Rabbit Monoclonal Antibody
Catalog #: AMRe12264



Summary

Production Name	HSPA12A (5U3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	HSPA12A FLJ13874 antibody;Heat shock 70 kDa protein 12A antibody;heat shock 70kD protein
Alternative Names	12A antibody;heat shock 70kDa protein 12A antibody;HS12A_HUMAN antibody;Hspa12a antibody;KIAA0417 antibody
Gene ID	259217.0
SwissProt ID	O43301.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	75kDa

Product Name: HSPA12A (5U3) Rabbit Monoclonal Antibody
Catalog #: AMRe12264




Background

HSPA12A is an ATP binding protein. Adapter protein for SORL1, but not SORT1. Delays SORL1 internalization and affects SORL1 subcellular localization.

Research Area

Image Data



Western blot analysis of HSPA12A expression in U-87 MG cell lysate.

Note

For research use only.