

Product Name: Liver Arginase (11G4) Rabbit Monoclonal Antibody
Catalog #: AMRe13341

Summary

Production Name	Liver Arginase (11G4) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	Arg1
Alternative Names	
Gene ID	29221.0
SwissProt ID	P07824.

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	35 kDa


Product Name: Liver Arginase (11G4) Rabbit Monoclonal Antibody
Catalog #: AMRe13341




Background

Research Area

Image Data



Western blot analysis of Liver Arginase in rat liver, PC-12 lysates using Liver Arginase antibody.



Western blot analysis of Liver Arginase in mouse liver lysates using Liver Arginase antibody.

Note

For research use only.