

Product Name: PCYT1A (4Z18) Rabbit Monoclonal Antibody
Catalog #: AMRe15871

Summary

Production Name	PCYT1A (4Z18) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	PCYT1A
Alternative Names	CCT A; CCT alph; CTPCT; PCYT1; Pcyt1a;
Gene ID	5130.0
SwissProt ID	P49585.A synthetic peptide of human PCYT1A

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	42kDa


Background

Product Name: PCYT1A (4Z18) Rabbit Monoclonal Antibody
Catalog #: AMRe15871

Controls phosphatidylcholine synthesis. Catalyzes the key rate-limiting step in the CDP-choline pathway for phosphatidylcholine biosynthesis.

Research Area

Image Data



Western blot detection of CT A in MCF7, HeLa cell lysates using CT A antibody (1:1000 diluted).

Note

For research use only.