

Product Name: Thioredoxin (6O1) Rabbit Monoclonal Antibody
Catalog #: AMRe18881



Summary

Production Name	Thioredoxin (6O1) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% New type preservative N and 0.05% BSA
Purification	Affinity purification

Immunogen

Gene Name	TXN
Alternative Names	
Gene ID	7295.0
SwissProt ID	P10599.

Application

Dilution Ratio	WB: 1/1000-1/10000
Molecular Weight	12 kDa


Product Name: Thioredoxin (601) Rabbit Monoclonal Antibody
Catalog #: AMRe18881




Background

Research Area

Image Data



Western blot analysis of Thioredoxin/TRX in C6, C2C12 lysates using Thioredoxin antibody.



Western blot analysis of Thioredoxin/TRX in HeLa, A549, HL-60, U2OS, U87-MG lysates using Thioredoxin/TRX antibody.

Note

For research use only.