

Product Name: TRIM56 (5Y14) Rabbit Monoclonal Antibody
Catalog #: AMRe19274



Summary

Production Name	TRIM56 (5Y14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	TRIM56
Alternative Names	
Gene ID	81844.0
SwissProt ID	Q9BRZ2.A synthetic peptide of human TRIM56

Application

Dilution Ratio	WB:1:10000
Molecular Weight	82kDa


Background

Product Name: TRIM56 (5Y14) Rabbit Monoclonal Antibody
Catalog #: AMRe19274



Research Area

Image Data



Western blot detection of TRIM56 in HeLa cell lysates using TRIM56 antibody(1:1000 diluted).

Note

For research use only.