

Product Name: TrkA (14C16) Rabbit Monoclonal Antibody
Catalog #: AMRe19289

Summary

Production Name	TrkA (14C16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	NTRK1
Alternative Names	NTRK1; MTC; TRK; TRK1; TRKA; Trk-A; p140-TrkA;
Gene ID	4914.0
SwissProt ID	P04629.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	87kDa


Background

The family of Trk receptor tyrosine kinases consists of TrkA, TrkB, and TrkC. While the sequence of these family members is highly conserved, they are activated by different neurotrophins: TrkA by NGF, TrkB by BDNF or NT4, and TrkC by NT3. Neurotrophin signaling through these receptors regulates a number of physiological processes, such as cell survival, proliferation, neural development, and axon and dendrite growth and patterning. Receptor tyrosine kinase involved in the development and the maturation of the central and peripheral nervous systems through regulation of proliferation, differentiation and survival of sympathetic and nervous neurons. High affinity receptor for NGF which is its primary ligand (PubMed: [1850821](http://www.uniprot.org/citations/1850821), PubMed: [1849459](http://www.uniprot.org/citations/1849459), PubMed: [1281417](http://www.uniprot.org/citations/1281417), PubMed: [8325889](http://www.uniprot.org/citations/8325889), PubMed: [15488758](http://www.uniprot.org/citations/15488758), PubMed: [22649032](http://www.uniprot.org/citations/22649032), PubMed: [17196528](http://www.uniprot.org/citations/17196528), PubMed: [27445338](http://www.uniprot.org/citations/27445338)). Can also bind and be activated by NTF3/neurotrophin- 3. However, NTF3 only supports axonal extension through NTRK1 but has no effect on neuron survival (By similarity). Upon dimeric NGF ligand- binding, undergoes homodimerization, autophosphorylation and activation (PubMed: [1281417](http://www.uniprot.org/citations/1281417)). Recruits, phosphorylates and/or activates several downstream effectors including SHC1, FRS2, SH2B1, SH2B2 and PLCG1 that regulate distinct overlapping signaling cascades driving cell survival and differentiation. Through SHC1 and FRS2 activates a GRB2-Ras-MAPK cascade that regulates cell differentiation and survival. Through PLCG1 controls NF-Kappa-B activation and the transcription of genes involved in cell survival. Through SHC1 and SH2B1 controls a Ras-PI3 kinase-AKT1 signaling cascade that is also regulating survival. In absence of ligand and activation, may promote cell death, making the survival of neurons dependent on trophic factors.

Research Area

Image Data

Product Name: TrkA (14C16) Rabbit Monoclonal Antibody
Catalog #: AMRe19289



Western blot analysis of TrkA expression in Human fetal brain tissue lysate.

Note

For research use only.