

Summary

Production Name	Flt3 Ligand Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	FLT3LG
Alternative Names	FLT3LG; Fms-related tyrosine kinase 3 ligand; Flt3 ligand; Flt3L; SL cytokine
Gene ID	2323
SwissProt ID	P49771

Application

Dilution Ratio	WB: 1/500-1/1000 ELISA: 1/10000
Molecular Weight	Calculated MW: 26 kDa; Observed MW: 26 kDa

Background

Stimulates the proliferation of early hematopoietic cells by activating FLT3.


Product Name: Flt3 Ligand Rabbit Polyclonal Antibody
Catalog #: APRab00627



Research Area

Cardiovascular

Image Data



Western blot analysis of Flt3 Ligand in HeLa lysates using Flt3 Ligand antibody.

Note

For research use only.