
Summary

Production Name	FGF11 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	FGF11
Alternative Names	FGF11; FHF3; Fibroblast growth factor 11; FGF-11; Fibroblast growth factor homologous factor 3; FHF-3
Gene ID	2256
SwissProt ID	Q92914

Application

Dilution Ratio	WB: 1/500-1/1000 ELISA: 1/10000
Molecular Weight	Calculated MW: 25 kDa; Observed MW: 25 kDa

Background

Product Name: FGF11 Rabbit Polyclonal Antibody
Catalog #: APRab00628




Probably involved in nervous system development and function.

Research Area

Cardiovascular

Image Data



Western blot analysis of FGF11 in NIH3T3, HeLa lysates using FGF11 antibody.

Note

For research use only.