

Summary

Production Name	PLD4 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	PLD4
Alternative Names	Phospholipase D4 (PLD 4; EC 3.1.4.4; Choline phosphatase 4; Phosphatidylcholine-hydrolyzing phospholipase D4)
Gene ID	122618
SwissProt ID	Q96BZ4

Application

Dilution Ratio	WB: 1/500-1/1000 ELISA: 1/10000
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 56 kDa

Background

Product Name: PLD4 Rabbit Polyclonal Antibody
Catalog #: APRab00812




5'->3' DNA exonuclease which digests single-stranded DNA (ssDNA).

Research Area

Cardiovascular

Image Data



Western blot analysis of PLD4 in various lysates using PLD4 antibody.

Note

For research use only.