
Summary

Production Name	Arkadia Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	RNF111
Alternative Names	RNF111; E3 ubiquitin-protein ligase Arkadia; RING finger protein 111
Gene ID	54778.0
SwissProt ID	Q6ZNA4.The antiserum was produced against synthesized peptide derived from human RNF111. AA range:901-950

Application


Dilution Ratio	WB 1:500 - 1:2000. ELISA: 1:40000.
Molecular Weight	110kD

Background


ring finger protein 111(RNF111) Homo sapiens The protein encoded by this gene is a nuclear RING-domain containing E3 ubiquitin ligase. This protein interacts with the transforming growth factor (TGF) -beta/NODAL signaling pathway by promoting the ubiquitination and proteosomal degradation of negative regulators, like SMAD proteins, and thereby enhances TGF-beta target-gene transcription. As a modulator of the nodal signaling cascade, this gene plays a critical role in the induction of mesoderm during embryonic development. Alternative splicing of this gene results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2012],function:Acts in the NODAL pathway of mesoderm patterning during embryonic development. Acts downstream AXIN1 as an E3 ubiquitin-protein ligase which promotes the ubiquitination of inhibitory SMADs such as SMAD7, induces their proteasomal degradation and thereby enhances the transcriptional activity of TGF-beta and BMP. Activates Smad3/Smad4-dependent transcription by triggering signal-induced SnoN degradation.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 RING-type zinc finger.,subcellular location:Upon TGF-beta treatment, translocates from nucleus to cytosol.,subunit:Interacts with SMAD6, SMAD7, AXIN1, AXIN2 and SKIL isoform SNON. Part of a complex containing RNF111, AXIN1 and SMAD7.,tissue specificity:Broadly expressed.,

Research Area

Image Data



Western blot analysis of RNF111 Antibody. The lane on the right is blocked with the RNF111 peptide.



Product Name: Arkadia Rabbit Polyclonal Antibody
Catalog #: APRab07140



Western blot analysis of the lysates from Jurkat cells using RNF111 antibody.

Note

For research use only.