


Summary

Name	WIF-1/Wnt inhibitory factor 1
Purity	Greater than 95% as determined by reducing SDS-PAGE
Endotoxin level	<1 EU/μg as determined by LAL test.
Construction	Recombinant Human Wnt Inhibitory Factor 1 is produced by our Mammalian expression system and the target gene encoding Gly29-Trp379 is expressed with a 6His tag at the C-terminus.
Accession #	AAH18037.1
Host	Human Cells
Species	Human
Predicted Molecular Mass	39.47 KDa
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM Glycine-HCl, 10% Sucrose, 5% Mannitol, 0.05% Tween80, pH3.5.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
Reconstitution	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE image



Background

Product Name: Recombinant Human WIF-1 (C-6His)
Catalog #: PHH1831



Alternative Names

Wnt Inhibitory Factor 1; WIF-1; WIF1

Background

Wnt Inhibitory Factor 1 (WIF1) is a secreted protein, which binds WNT proteins and inhibits their activities. WNT proteins are extracellular signaling molecules involved in the control of embryonic development. WIF1 contains a WNT inhibitory factor (WIF) domain and 5 epidermal growth factor (EGF)-like domains. is found to be present in fish, amphibia and mammals. WIF1 is a recurrent target in human salivary gland oncogenesis. WIF1 may be involved in mesoderm segmentation. WIF1 is a tumor suppressor, specifically in nonfunctioning pituitary tumors.

Note

For Research Use Only , Not for Diagnostic Use.